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**10% NEUTRAL BUFFERED FORMALIN**

Ready-to-use solution containing 4% formaldehyde (weight/volume)

Single-use IN VITRO DIAGNOSTIC Medical Device

Basic UDI-ID: 37015164FORMOL4KJ

Manufacturer: Laurypath, 52 route de Brignais - 69630 Chaponost. FRANCE Tel: +33(0)4 78 34 21 67
www.laurypath.fr. Email: contact@laurypath.fr.

ISO 13485 certified quality management system.

Class A in vitro diagnostic device according to regulation 2017/746

1. INTENDED USE

Intended to be used as a liquid conservation and transport medium for biopsies and surgical specimens in diagnostic procedures.

REFERENCE	DESCRIPTION
STANDARD PREFILL POTS	
40877-02/10E	24 prefill pots containing 10 ml
40877-02/10EB	24 prefill pots containing 10 ml
40877-02/10ER	24 prefill pots containing 10 ml
40877-02/10EV	24 prefill pots containing 10 ml
40877-04/20E	24 prefill pots containing 20 ml
40877-04/20EB	24 prefill pots containing 20 ml
40877-04/20ER	24 prefill pots containing 20 ml
40877-04/20EV	24 prefill pots containing 20 ml
40877-06N30ER	35 prefill pots containing 30 ml
40877-06N30E	35 prefill pots containing 30 ml
40877-06N30EB	35 prefill pots containing 30 ml
40877-06N30EV	35 prefill pots containing 30 ml
40877-125/60E	24 prefill pots containing 60 ml
40877-125/60EB	24 prefill pots containing 60 ml
40877-125/60ER	24 prefill pots containing 60 ml
40877-125/60EV	24 prefill pots containing 60 ml
40877-180/80E	24 prefill pots containing 80 ml
40877-180/80EB	24 prefill pots containing 80 ml
40877-180/80ER	24 prefill pots containing 80 ml
40877-180/80EV	24 prefill pots containing 80 ml
CUSTOMISED SPECIFIC PREFILL POTS	
F-02/10EX-XXX	24 prefill pots containing 10 ml
F-04/20EX-XXX	24 prefill pots containing 20 ml
F-06N30EX-XXX	35 prefill pots containing 30 ml

REFERENCE	DESCRIPTION
CANISTERS	
40877-29 - XXX	1-L Canister
40877-36 - XXX	5-L Canister
40877-41 - XXX	10-L Canister
40877-46 - XXX	25-L Canister
STANDARD PREFILL CONTAINERS	
40877-26/250E - XXX	Prefill Container containing 250 ml
40877-29/500E - XXX	Prefill Container containing 500 ml
SPECIFIC PREFILL CONTAINERS	
F-26/250E-XXX	Prefill Container containing 250 ml
F-29/500E-XXX	Prefill Container containing 500 ml
STANDARD PREFILL CONTAINERS	
40877-33/1500E - XXX	Prefill Container containing 1500 ml
40877-36/2500E - XXX	Prefill Container containing 2500 ml
CUSTOMISED PREFILL CONTAINERS	
F-33/1500E-XXX	Prefill Container containing 1500 ml
F-36/2500E-XXX	Prefill Container containing 2500 ml

F-125/60EX-XXX	24 prefill pots containing 60 ml
F-180/80EX-XXX	24 prefill pots containing 80 ml

XXX: 3 optional letters in the case of a customised product label.

2. PRINCIPLE OF THE METHOD

10% Neutral Buffered Formalin is used as a fixative in histology and cytology. The presence of phosphate buffer enables a pH of 6.8 – 7.2 to be maintained, avoiding the formation of formalin pigments. Methanol is used to limit the polymerisation of formaldehyde.

3. COMPOSITION

10% neutral buffered formalin is expressed by volume of concentrated formalin on the total volume of solution. It is also called buffered 4% formaldehyde if expressed by weight of formaldehyde on the volume of solution. Its pH is between 6.8 and 7.2.

Formaldehyde CAS No. 50-00-0

Methanol CAS No. 67-56-1

Monopotassium phosphate CAS No. 7778-77-0

Disodium phosphate CAS No. 7558-79-4

4. TREATMENT BEFORE USE

Ready-to-use solution.

5. QUALITY CONTROL

All production batches are released after verification of their conformity through quality control inspections and analyses. Batch numbers of the products and the raw materials that they contain are recorded using validated computerised systems that guarantee traceability.

6. INSTRUCTIONS FOR USE

The product is intended exclusively for professional use and must be used by qualified and experienced personnel, prefill pots are designed for single use only.


1. Rinse samples covered with organic matter (blood, bile, faeces, etc.), as this can slow down fixation and compromise microscopic interpretations.
2. Fully immerse the sample in formalin. Sample/fixative volumetric ratio 1:10. Maximum sample thickness: 1 cm.
3. Comply with the laboratory protocols. The fixation time depends on the thickness and shape of the samples. The average penetration speed of formalin is 1 mm per hour.
4. Do not reuse solutions that have been used.

NOTE: If used in enclosed tissue processors it is advisable to follow the manufacturers instructions regarding flushing and housekeeping procedures.

For containers and pots pre-filled with formalin, ensure they are closed properly before transporting them to avoid any leakage of formalin or loss of biopsies.

To close the pots correctly, follow these instructions:

- 1 - Place the cap (lid) on the container and turn once anti-clockwise. You will hear a noise which corresponds to the correct positioning of the cap on the container.

2 - Close the container by turning the cap clockwise and fully tightening the cap.

To close the container correctly, follow the instructions below:

7. STORAGE

Store the reagent in original packaging at a temperature between 0° and 40°C. Do not freeze and avoid exposing to direct sunlight. When correctly stored in its undamaged packaging, the product is stable until the expiry date indicated on the label. The Prefill Containers must be properly closed.

8. DISPOSAL INSTRUCTIONS

Dispose of the product in an appropriate hazardous product reprocessing channel.

9. WARNINGS AND PRECAUTIONS FOR USE

Do not use this product after its expiration date dictated on the product label.

The product is ready to use. However, if it is transferred to another Prefill Container, report its reference, batch number, expiration date, hazard statements and danger pictograms.

For any pressurised (pneumatic) transport use, the compatibility of the devices with these systems must be checked beforehand.

Before introducing the sample, check that the liquid is not cloudy or does not have any traces of a deposit, if it is do not use the product and contact the manufacturer.

10. SAFETY

Read the label and the safety data sheet before use.

The safety data sheet is available on request or on our website www.laurypath.fr.

Harmful by inhalation, in contact with skin and if swallowed.

To be used exclusively in ventilated conditions (suction, hood, etc.).

Wear personal protective equipment and do not inhale fumes.

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user/or the patient resides.

Human samples are potentially biohazardous. Good Laboratory Practice must be used when handling tissue samples to prevent cross contamination and infection. The user should complete a risk assessment to determine any potential hazards related to tissue handling.

11. SYMBOL KEY

Batch No.

Storage temperature limits

Manufacturer

Product reference

In vitro diagnostic medical device

Expiration date

Single-use device

Read the instructions for use carefully

Unique identification of the medical device

12. MD UNIQUE IDENTIFICATION

(01) 0 3701516 40067 0
(17) 240529
(10) 227952

(01): IUD – ID: Unique identification of the medical device, specific to each reference

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